

# Vasco® Surgical Powder-free

## STERILE SURGICAL AND PROTECTIVE GLOVES | DATA SHEET



**B. Braun Melsungen AG confirms that  
Vasco® Surgical Powder-free gloves comply with the following standards,  
directives and regulations:**

### EC CERTIFICATES AND APPLIED STANDARDS

Medical Device Class IIa CE 0123 (TÜV Süd, DE), according to MDD 93/42/EEC

EN 455 1-4, ISO 10282, ISO 10993, ISO 11137

ASTM D3577, ASTM D5712, ASTM D6978

Personal Protective Equipment Category III according to Personal Protective Equipment Regulation (PPER) EU 2016/425

EN 421, EN 420, EN 374, ISO 16523, ISO 16604, ASTM F1671

ISO 9001, ISO 13485

### QUALITY CERTIFICATES

### PERSONAL PROTECTIVE EQUIPMENT

Information and Declaration of Conformity according to PPER (EU) 2016/425:



[www.bbraun.com/gloves-declarations-of-conformity](http://www.bbraun.com/gloves-declarations-of-conformity)

B. Braun Melsungen AG

A handwritten signature in blue ink, belonging to Dr. Hans-Ulrich Gaudin, positioned above his printed name and title.

Dr. Hans-Ulrich Gaudin  
Head of Global Regulatory Affairs OPM Germany

# Vasco® Surgical Powder-free

## STERILE SURGICAL AND PROTECTIVE GLOVES | REGULATORY INFORMATION

### MEDICAL DEVICE INFORMATION


MDD 93/42/EEC (CLASS IIa), EN 455



### PERSONAL PROTECTIVE EQUIPMENT INFORMATION

Tested in accordance with:

ISO 374-1:2016/Type B



KPT

CE 2777

PPE Regulation (EU) 2016/425 (Cat. III);  
EN 420:2003+A1:2009

Code letter	Test chemical	EN 374-1:2016 Permeation level	EN 374-4:2013 Mean degradation
K	Sodium hydroxide 40 %	Level 6	-11,0 %
P	Hydrogen peroxide 30 %	Level 6	10,1 %
T	Formaldehyde 37 %	Level 6	2,6 %

Tested acc. to EN 16523-1:2015

Performance levels acc. EN 374-1:2016 +A1:2018	1	2	3	4	5	6
Measured breakthrough times (mins)	> 10	> 30	> 60	> 120	> 240	> 480

Degradation levels indicate the change in puncture resistance of the gloves after exposure to the challenge chemical. NOTE: Where the test specimens gave an increased puncture force after chemical exposure, the result is reported as a negative degradation.

AQL 0.65

Resistance to bacteria and fungi	pass
Resistance to virus	pass

ISO 374-5:2016



VIRUS

EN 421:2010



Protection against particulate radioactive contamination.

This information does not reflect the actual duration of protection in the workplace and the differentiation between mixtures and pure chemicals. The chemical and penetration resistance has been assessed under laboratory conditions from samples taken from the palm only and relates only to the chemical tested. It can be different if the chemical is used in a mixture. It is recommended to check that the gloves are suitable for the intended use because the conditions at the workplace may differ from the type test depending on temperature, abrasion and degradation. When used, protective gloves may provide less resistance to the dangerous chemical due to changes in physical properties. Movements, snagging, rubbing, degradation caused by the chemical contact etc. may reduce the actual use time significantly. For corrosive chemicals, degradation can be the most important factor to consider in selection of chemical resistant gloves. Before usage, inspect the gloves for any defect or imperfections.

# Vasco® Surgical Powder-free

## STERILE SURGICAL AND PROTECTIVE GLOVES | TECHNICAL DATA



SIZE	REF	GLOVE DIMENSIONS (EN 455)	
		Width of palm	Total length
5.5	6081100	70 ± 4 mm	≥ 290 mm
6	6081101	77 ± 5 mm	
6.5	6081111	83 ± 5 mm	
7	6081121	89 ± 5 mm	
7.5	6081131	95 ± 5 mm	
8	6081141	102 ± 6 mm	
8.5	6081151	108 ± 6 mm	
9	6081161	114 ± 6 mm	

### PHYSICAL PROPERTIES

		Min. specification	Typical value
Wall thickness	Palm	0.21 mm	0.21 mm
	Cuff	0.17 mm	0.17 mm
Force at break (acc. to EN 455)	During shelf life	9 N	15 N before ageing 13 N after ageing
Elongation at break (acc. to ASTM D 3577)	Before ageing	750 %	891 %
	After ageing	560 %	869 %
Tensile strength (acc. to ASTM D 3577)	Before ageing	24 MPa	29 MPa
	After ageing	18 MPa	26 MPa

### GLOVE DESIGN

Colour	natural white
Shape	fully anatomical shape with curved fingers
Cuff	rolled rim
Surface finish	textured palm and finger area, silicone treated
Inner glove surface	polymer coated, powder-free

### GLOVE MATERIAL

Natural rubber latex	protein content < 50 µg/g
Latex allergy risk	containing natural rubber latex which may cause allergic reactions including anaphylactic reactions

### ACCELERATORS

Zn-dithiocarbamate	
Free of thiurames, thioureas and thiazoles - including mercaptobenzothiazole MBT	

### STERILIZATION

Electron beam radiation
-------------------------

### LOGISTIC INFORMATION

Peel pouch	1 pair	270 x 140 mm (L x W)
Dispenser pack	50 pairs	270 x 140 x 215 mm (L x W X H)
Transportation carton	4 dispenser packs	548 x 147 x 439 mm (L x W X H)
Shelf life	5 years	
Storage conditions	store at room temperature, protect from dust, humidity, sun light and ozone	

# Vasco® Surgical Powder-free

## STERILE SURGICAL AND PROTECTIVE GLOVES | BARRIER PROPERTIES – CYTOSTATIC DRUGS



### CLASSIFICATION

- Not suitable
- Suitable if changed before permeation breakthrough
- Suitable for prolonged use

Tested by ARDL, USA in accordance with

**ASTM D 6978:** Standard Practice for Assessment of Resistance of Medical Gloves to Permeation by Chemotherapy Drugs. Minimum detection rate 0,01 µg/cm<sup>2</sup>/min

CHEMOTHERAPY DRUG	mg/ml	CAS registry no.	MIN BREAKTHROUGH DETECTION TIME	
Carmustine	3.3	154-93-8	<span style="color: teal;">■</span>	> 240 min
Cisplatin	1.0	15663-27-1	<span style="color: teal;">■</span>	> 240 min
Cyclophosphamide monohydrate	20.0	6055-19-2	<span style="color: teal;">■</span>	> 240 min
Doxorubicin hydrochloride	2.0	25316-40-9	<span style="color: teal;">■</span>	> 240 min
Etoposide	20.0	33419-42-0	<span style="color: teal;">■</span>	> 240 min
Fluorouracil	50.0	51-21-8	<span style="color: teal;">■</span>	> 240 min
Mitoxantrone	2.0	65271-80-9	<span style="color: teal;">■</span>	> 240 min
Paclitaxel (Taxol)	6.0	33069-62-4	<span style="color: teal;">■</span>	> 240 min
Thio-Tepa	10.0	52-24-4	<span style="color: teal;">■</span>	> 240 min