Ringer B. Braun

BG	Инструкции за употреба	2
CZ	informace o použití	3
DE	Gebrauchsanweisung	4
DK	Brugsanvisning	5
EE	Kasutusteave	6
ES	Instrucciones de uso	7
FI	Tietoa käyttäjälle	8
FR	Mode d'emploi	9
GB	Information for use	10
GR	Πληροφορίες σχετικά με τη χρήση	11
HU	Használati utasítás	12
IT	Foglietto illustrativo	13
LT	informacija apie naudojimą	14
LV	informācija lietošanai	15
NL	Gebruiksaanwijzing	16
NO	Bruksanvisning	17
PL	Sposób użycia	18
PT	Informações para a utilização	19
RO	Informații pentru utilizare	20
SE	Bruksanvisning	21
SI	Navodilo za uporabo	22
SK	Informácia o použití	23

schwarz

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Lätus - - -

Irrigation Solution___191_GIF
191/12609812/0413
GIF (A5-Booklet, 22-spr.)
Standort Crissier

B | BRAUN

Ringer B. Braun

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GB – Information for use

Ringer B. Braun

Irrigation solution,
sterile and endotoxin-free

Composition

1000 ml solution contain:

Sodium chloride	8.60 g
Potassium chloride	0.30 g
Calcium chloride 2H ₂ O	0.33 g
in Water for Injection	
Electrolytes:	mmol/l
Sodium	147.20
Potassium	4.02
Calcium	2.24
Chloride	155.70

Indications

For irrigation of wounds and burns. For irrigation and cleansing during surgical intervention, during endoscopic examinations of body cavities, intra-operative arthroscopy with mechanical instrumentarium and for postoperative irrigation.

Dosage

The amount of fluid depends on the clinical condition.

Overdosage

during wound irrigation and intra-operative application can lead to absorption of irrigation fluid, that can lead to hypervolaemia in the form of hypotonic hyperhydration (hyponatraemia). Headache, nausea, restlessness and disorientation are symptoms of hyperhydration. In severe cases coma-like states can occur.

Countermeasures in cases of overdosage

Discontinuation of irrigation, correction of the water and electrolyte balance depending on the particular clinical situation.

Warnings and precautions

Don't use for infusion. Only use the solution if the seal is undamaged and the solution is clear. Opened containers must not be saved, quantities left over must be discarded. Re-use of single-use devices creates a potential risk of patient or user. It may lead to contamination and/or impairment of functional capability. Contamination and/or limited functionality of the device may lead to injury, illness or death of the patient. Store out of reach of children. When there is a risk of absorption of larger amounts

of the solution, as may be the case when using it in large body cavities or large wound surfaces, it should not be used in patients in states of hyperhydration.

Under those conditions it should be used with caution in patients with

- Hypertonic dehydration,
- Hyperkalaemia
- Hyponatraemia,
- Hyperchloraemia,
- Renal insufficiency with a tendency to hyperkalaemia,
- Disorders where restriction of sodium intake is indicated, such as cardiac insufficiency, generalised oedema, pulmonary oedema, hypertension, eclampsia, severe renal insufficiency.

Monitoring should include checks of the serum electrolytes and the water balance. Ringer B. Braun must not be used for intrathecal irrigation procedures.

Adverse side effects

None are known.

Expiry date

The solution must not be used after the expiry date given on the label.

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B. Braun Melsungen AG