

KRUUSE Lite-Absorb

Absorbent dressing with light absorbency properties for use in mildly exudating wounds.

Ideal for use with KRUUSE Hydrogel or KRUUSE Manuka Honey which both contribute to keeping the wound moist.

Application

- Use as a primary layer dressing to cover wounds and absorb wound exudate or as a secondary layer to manage a variety of wound conditions
- Easily secure the dressing with a gauze bandage or tape

Improved healing

- The dressing is made of viscose and polyester and has a perforated PE film on both sides reducing the risk of attachment to the wound bed
- The tightly woven layers ensure removal of the bandage in one piece, without leaving lint in the wound

More comfort

 The cushioned inner layer protects the wound, and the non-adherent PE film (polyethylene) ensures decreased pain during dressing change and removal

Cat. No 160315 160316 160317	Description KRUUSE Lite-Absorb KRUUSE Lite-Absorb KRUUSE Lite-Absorb	Size, Metric 5 x 5 cm 10 x 10 cm 10 x 20 cm	Size, Imperial 2" x 2" 4" x 4" 4" x 8"	Package 25/pk 25/pk 25/pk	
		Absorbent vi	on-adherent PE film scose and polyester on-adherent PE film		
		ALCON TO SERVICE STATE OF THE			
					Scan QR code to watch the video