

BUSTER Green Ocean Collar

Rethinking pet protection

We transform discarded fishing nets into E-collars so that you can make the sustainable choice that protects pets and the world we live in. This is our initiative to rethink, recycle and provide environmentally friendly products for your practice.

*BUSTER protective solutions
just got more protective:
For pets
For the oceans
For marine life worldwide*




Recycling: With the BUSTER Green Ocean Collar, we deliver superior pet protection and empower you to make the natural choice for you practice and the wider world.

Why go green?

To maximise old materials

Recycling old fishing nets is a new focus area in the booming plastics recycling industry. Without it, this plastic would be sent to landfill, burnt, or lost at sea - destroying habitats and endangering marine animals.

The BUSTER Green Ocean Collar is made from old fishing nets that are collected and recycled responsibly.



RETHINK
AND
RECYCLE



MAXIMISE
OLD
MATERIALS



EMPOWER
GREEN
CHOICES

From disregarded to liberated

Discarded fishing nets in the ocean threaten marine animals, seabirds, and the entire marine environment.

Imagine the weight of over 3,000 fully loaded jumbo jets ready for takeoff. This is how much plastic waste from the maritime industry ends up in our environment every year.

Up to 13 million tons of plastic waste
end up in our environment and oceans every year.

Around 10% of this plastic waste
stems from discarded fishing nets from the maritime industry.

If this continues, there could be more plastic
than fish in the oceans by 2050.



The facts

- Fishing activities are intensifying around the world
- Discarded fishing nets are increasing significantly
- Plastic pollution is a huge problem affecting us all
- We can all act to do something about it

Sources: <https://bit.ly/3zpdduH>

Same reliable standard

The BUSTER range has become known as the benchmark in pet protection and the BUSTER Green Ocean Collar is no different. It is equipped with the easy-to-apply “Clic” closing system with outer straps for secure fastening to the dog’s neck collar. The opaqueness of the collar provides a den-like and stressless environment for the dog.

Perfect for

- Postoperative care
- Wound care
- Spay/neuter recovery
- Hot spots and skin diseases

Provides

- Fast recovery
- Full body coverage
- Optimal security

Ensures

- Stressless handling
- User-friendliness
- Pet owner compliance



BUSTER Green Ocean Collar

Cat. No	Size	Length in cm / Inch	Opening in cm / Inch	Neck-hole in cm / Inch
273342	7.5 cm	7.5 cm / 3"	15-17 cm / 6"-7"	7-8 cm / 3"
273343	10 cm	10 cm / 4"	17-19 cm / 7"	7-8 cm / 3"
273344	12.5 cm	12.5 cm / 5"	19-23 cm / 7"-9"	9-10.5 cm / 4"
273345	15 cm	15 cm / 6"	24-28 cm / 9"-11"	10-12 cm / 4"-5"
273346	20 cm	20 cm / 8"	32-35 cm / 13"-14"	10-13 cm / 4"-5"
273347	25 cm	25 cm / 10"	36-43 cm / 14"-17"	12-15 cm / 5"-6"
273348	30 cm	30 cm / 12"	42-51 cm / 17"-20"	15-18 cm / 6"-7"
273349	Starter Pack of sizes 7.5 - 30 cm			